

This medical device is intended for patient positioning.

This is not a restraint system.

*Before use, it is important to check the status of the device: stitching, state of the buckles, slings, fabric.
Any damage will either be repaired by NAUSICAA Medical, or will result in a disposal.*

ECO BELT

FOR ALL CHAIRS

Abdominal ECO Belt: CECO-A

Pelvic ECO Belt: CECO-P

Bust ECO Belt: CECO-B

Integral ECO Belt: CECO-I



CECO-A



CECO-P



CECO-B



CECO-I



**THIS MEDICAL DEVICE MUST BE
IMPLEMENTED BY A TRAINED
HEALTH PROFESSIONAL**

OPTIONS

Pelvic
Modular
Element



HOME CARE RANGE



CECO-A	ABDOMINAL ECO BELT	One Size Fits	<input type="radio"/>
CECO-P	PELVIC ECO BELT	One Size Fits	<input type="radio"/>
CECO-B	BUST ECO BELT	One Size Fits	<input type="radio"/>
CECO-I	INTEGRAL ECO BELT	One Size Fits	<input type="radio"/>

Lot

PACKAGING : 1 PIECE

Textile materials in
accordance with the legislation
REACH
on stable substances

12 Allée du Piot - ZA Pôle Actif
30660 Gallargues le Montueux
FRANCE

CE
v.04.03.2017

IMPORTANT RECOMMENDATIONS :

In order to get a maximum efficiency of these sort of product, it is necessary:

- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient

These products must be in direct contact with a wounded skin.

ECO BELT

COMPOSITION

HEMOCARE RANGE

Jersey polyester doublé de mousse polyuréthane

Biais 91 % polyamide - 9 % élasthanne

Sangle 100 % polyester

Boucles nylon

Épaulières PVC

INSTRUCTION FOR USE

The ECO Belt is declined under 4 modular versions :

1. the Abdominal ECO Belt.

-Position the belt at the patient's waist level.

- Pass the sling behind the backrest of the chair

- Close the nylon buckle, then adjust the sling length by avoiding compress the patient's abdominal

2. the Pelvic ECO Belt is composed of the Abdominal ECO Belt to which is added a pelvic modular element.

-Position the pelvic maintain on the chair, fix the sling 25 mm at the backrest of the chair

-Seat the patient on the chair

-Position the slings behind the behind the headrest of the chair. Close the nylon buckle, then adjust the sling length by avoiding compress the patient's abdominal.

3. the Bust ECO Belt.

-Position the belt at the patient's waist level.

- Pass the sling behind the backrest of the chair

-Close the nylon buckle, then adjust the sling length by avoiding compress the patient's abdominal

-Set the shoulder harness by positioning correctly the shoulder pad on the patient's shoulders, put the buckle behind the headrest. Pass a shoulder pad sling under the maintaining sling, connect the buckle and adjust its length

4. the Integral ECO Belt is composed of the Bust ECO Belt to which is added a pelvic modular element.

For the ECO Belt Integral settlement, proceed firstly like the Pelvic ECO Belt (2) then like the paragraph d) concerning Bust ECO belt (3).

v 04.03/2017



**CE DISPOSITIF MÉDICAL DOIT ÊTRE MIS EN ŒUVRE OU FAIRE
L'OBJET D'UNE FORMATION PAR UN PROFESSIONNEL DE SANTÉ**