

STAND-UP LIFTS SLINGS RANGE USER GUIDE

Size Chart

SIZES	TU (one size fits)	S	M	L	XL
BUST (in cm)	/	70/90	90/115	115/140	140/165
WAIST (in cm)	/	70/90	90/115	115/140	140/165



Choosing the correct sling size is essential in order to ensure safety during the transfer.

A wrong sized sling generates discomfort and increases the risk of injury for the patient.

A size colour-code of the edges has been put in place in order to quickly identify the size of each sling:

- to the size S corresponds the Yellow colour
- to the size M corresponds the Red colour
- to the size L corresponds the Green colour
- to the size XL corresponds the Blue colour
- to the size TU (one size fits) corresponds the Black colour

Compatibility Chart

	WAYUP 5 (page 14)	WAYUP 5 with opening legs (page 15) & WAYUP XL (page 16)	EASYLEV 5 (page 20)	EASYLEV 5 with opening legs (page 21) & EASYLEV XL (page 22)	BLUE WAYUP (page 26)	BLUE WAYUP with opening legs (page 27)	BLUE EASYLEV (page 32)	BLUE EASYLEV with opening legs (page 33)
ECO SLING (page 38)	★	★	★	★	★★★★	★★★★	★★★★	★★★★
SPECIFIC SLING (page 38)	★★★★	★★★★	★★	★★	★★★★	★★★★	★★	★★
STANDARD SLING (page 39)	★★	★★	★★★★	★★★★	★★	★★	★★★★	★★★★

In order to select the size of a sling, the caregiver must match the chest (Standard and Eco Slings) or waist (Specific Sling) with the dimensions mentioned in the table beside.

IMPORTANT

An assessment of the patient's capacity is required before using a stand-up lift; should be considered, depending on the situation, the patient's physical handicap, pathology, and morphology.