

MODULUP PRODUCT SHEET





Description :

- Frame made in steel
- Laser cutting and bending sheet metal, tubes and profiles
- Epoxy baked painting
- Ball bearing castors
- Weight capacity : 150 kg
- Lifespan : 5 years

- Compliance with the following standards :
 - Medical products (93/42 EEC amended by 2007/47/EEC)
- Applicable standards :
 - BS EN ISO 10535 : 2007 (Hoists for transfer of disabled persons-Requirements and test methods)

Dimensionals Characteristics :



Dimensions (in mm)		
1	Highest position	/
2	Lowest position	/
c	CSP Maximum height	1030
d	CSP Minimum height	540
f	Maximum intern width	580
g	Minimum intern width at maximal position	390
h	Minimum intern width	340
i	Garde au sol	34
j	Base height	80
k	Minimum distance between the CSP and the wall at maximum height	480
m	Minimum distance between the CSP and the wall at minimum height	290
n	Turning diameter	860