

-  300 kg
-  46 kg
-  660/1820 mm
-  113 mm
-  700/1050 mm




SCAN ME
to learn more

FLYER XL

electrical opening legs

It is a high-end device designed to answer the requirement that bariatric quality treatment must be carried out in specific devices as it is powerful due to its lifting actuator (10 000 N) and convenient with its electrical opening legs (up to 105 cm) as standard.

WATCH
the product video on
our YOUTUBE channel




VIDEO



SCAN
to watch





Ergonomically designed 4-point spreader bar for an optimized weight distribution



10 000 N actuator for an increased maximum weight capacity (300 kg)



Premium control box for an enhanced use



Strong casters with double ball bearing

OPTIONS



Lithium Battery



Weighing System

REFERENCES
with standard battery: FLYER-XL-E-PF
with lithium battery: FLYER-XL-E-LI-PF