



SCAN ME
to learn more

HOLST ARM SLING

It is designed to prevent or to completely reduce the glenohumeral sub-luxation (source of pains) while allowing a daily setting up by the patients themselves, without changing the initial adjustment.

The elasticated strap exerts a vertical traction on the humerus by applying the force of the traction to the elbow and to the homolateral shoulder.



REFERENCES

LEFT SIDE: **EC-HSTG-PF**
RIGHT SIDE: **EC-HSTD-PF**



SCAN ME
to learn more

ARM SLINGS

Available in 2 models (with or without a dorsal strap), they are used in case of fracture, tendonitis or sprain to promote healing of the upper limb or shoulder.

Esthetical and comfortable, they are also equipped with a thumb rest.

The addition of a dorsal strap keeps the sling against the user's body (dorsal strap model only).



Arm sling



Arm sling with dorsal strap

REFERENCES

ARM SLING
Size 1: **ECB-01-PF**
Size 2: **ECB-02-PF**
Size 3: **ECB-03-PF**

ARM SLING WITH DORSAL STRAP
Size 1: **ECBSM1-PF**
Size 2: **ECBSM2-PF**
Size 3: **ECBSM3-PF**