


# PRESENTATION



Our Patient Positioning range participates in the search for a good posture of patients while preventing falls, is CE certified and is in compliance with REACH legislation.

The purpose of positioning is to provide a stable and comfortable support for patients confined to a chair or a bed, as its design makes it easy for everyday use and the choice of fabrics makes it comfortable, while being safe.

Our range is developed in partnership with health professionals (occupational therapists and physiotherapists) in order to ensure the best possible support for the patient as well as a quick and easy installation for the caregiver.

The following products are available in either Home Care Range, Hospital Range or both.

The Home Care Range products (washable at 30°C) are made of a foam lamination composed of bengaline/cotton polyester on one face and of jersey/polyester fabric on the other face.

The Hospital Range products (washable at 70°C) are manufactured with a lamination of foam, and on both faces, of a waterproof and permeable jersey/PU coating.

## CE CERTIFICATION

All our products are in compliance with EU requirements arising to the product's manufacturer.

## WASHING

The Home Care Range products are washable at 30°C.  
The Hospital Range products are washable at 70°C.

## REACH COMPLIANCE

Our slings are in compliance with the REACH regulation which aims to progressively remove the most dangerous chemical substances in the European Union.

## COMPATIBILITY

The NAUSI-SURE system is a secure restraining device with a magnetic closure that prevents the patient, or any unauthorized person to withdraw the product used.  
A magnetic key is available to the caregiver in order to quickly unlock the device to release the patient.



HOME CARE RANGE



HOSPITAL RANGE