



Prior to each use, it is important to check the status of the device: seams, state of the straps and of the fabric. Any damage will either be repaired by NAUSICAA Medical, or will result in a disposal. For hygienic reasons and respect for health protocol, textile products have been opened will not be returned nor exchanged.

WALKING AID & LIFTING BELT



THIS MEDICAL DEVICE
MUST BE IMPLEMENTED
BY A TRAINED HEALTH
PROFESSIONAL



Equipped with a self-centering
50 mm buckle

Composition :
Combination of polyurethan foam
and jersey/cotton polyester fabric
Jersey/PU coating
Polyester strap
100% polyester fabric


Uses :

As a walking aid, it is used on the patient for a safe handling.

As a lifting device, it is used around the patient's waist to provide a gripping point to the caregiver.

TO BE USED ONLY FOR THE STATED PURPOSE:

- QUESTION-ON:
- THIS BELT CAN BE USED ON THE PATIENT FOR A MAINTENANCE SAFELY.
- LINKAGE:
- IT IS USED AROUND THE PATIENT'S WAIST TO PROVIDE A GRIPPING POINT TO THE CAREGIVER.



CADM2-PF