



This medical device is intended for patient positioning.

This is not a restraint system.

Before use, it is important to check the status of the device: stitching, state of the buckles, slings, fabric. Any damage will either be repaired by NAUSICAA Medical, or will result in a disposal.

ECO SUPPORT BELT

FOR ALL CHAIRS



THIS MEDICAL DEVICE MUST BE
IMPLEMENTED BY A TRAINED
HEALTH PROFESSIONAL




TO BE USED ONLY FOR THE
INDICATED USE :

**TO PREVENT FROM
FALLING**

Maximal strap length : + 164 cm
Belt height : + 11 cm
Buckle : 50 mm

Simple and one size fits all, it is designed for different chair uses.

- Maximum useful area: 164 cm
- Without size constraints



CCMED2HP-PF

Important Recommendations :

In order to get a maximum efficiency of these sort of product, it is necessary :

- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient

These products must not be in direct contact with a wounded skin.



ECO SUPPORT BELT

HOSPITAL RANGE

Combination of polyurethan foam
and on both faces
Jersey/PU coating
100 % polyester bias
100 % polyester strap
Nylon buckle

INSTRUCTIONS FOR USE

TO BE USED ONLY FOR THE INDICATED USE :
TO PREVENT FROM FALLING

- Place the belt around the patient's abdomen (at hip level) and close the loop on the back of the chair.

v 05.12/2020



THIS MEDICAL DEVICE MUST BE SET UP BY A PROFESSIONAL
OR THE USER MUST BE EDUCATED TO SET IT UP.