

ISO 13485

BUREAU VERITAS
Certification



NAUSICAA Médical

www.nausicaa-medical.eu

MADE IN FRANCE
IN OVER WORKSHOP
SINCE 1992



Prior to each use, it is important to check the status of the device: seams, state of the straps and of the fabric. Any damage will either be repaired by NAUSICAA Medical, or will result in a disposal. For hygienic reasons and respect for health protocol, textile products have been opened will not be returned nor exchanged.

LIFTING BELT

TO BE USED ONLY FOR THE
STATED PURPOSE :
AID TO LIFT A PERSON
SITTING




THIS MEDICAL DEVICE
MUST BE IMPLEMENTED
BY A TRAINED HEALTH
PROFESSIONAL



Jersey/PU coating
PE foam 100 kg/m³
Polyester strap
100% polyester fabric

Description

The Lifting Belt is used to raise a sitting person and preserves the caregiver of awkward postures during manual transfers. Its inside fabric is in anti-slipping PU coating.



CREL-PF

v 05.12/2020

