




*Prior to each use, it is important to check the status of the device:
seams, state of the straps and of the fabric.
Any damage will either be repaired by NAUSICAA Medical, or will result in a disposal.*

SLIDE SHEETS



**THIS MEDICAL DEVICE MUST BE
IMPLEMENTED BY A TRAINED
HEALTH PROFESSIONAL**



DGL-PF



SLIDE SHEET

HOSPITAL RANGE
Tissu 100% polyamide :
- Enduction silicone (dessus)
- Enduction polyuréthane (dessous)

DESCRIPTION

- For the handling of a patient, the Slide Sheet is a device that facilitates lateral and vertical translation as well as the raising of a lying patient in a bed.
- The sliding effect facilitates the work of carer, promotes the prevention of pressure ulcers and improves patient comfort.
- The Slide Sheet is composed a of synthetic fabric polyamide tube in polyurethane induction, the two layers of fabric in contact slide on one another.
- Silicone Treatment on Patient Side.

USES



v 05.12/2020



THIS MEDICAL DEVICE MUST BE SET UP BY A PROFESSIONAL
OR THE USER MUST BE EDUCATED TO SET IT UP.