


ECO SLING




Specially designed to be used with the BLUE WAYUP range (page 24) and the BLUE EASYLEV range (page 30), the Eco Sling is available in one size fits all.

Placed under the armpits or on the lumbar region, its padding under the arms distributes traction on the patient and avoids painful pressure during a transfer.



TU	A	107
	B	29

	REFERENCE
Jersey	SXPECO-PF-TU



SPECIFIC SLING




Manufactured to be used with the WAYUP 5 range (page 12), the Specific Sling is available in 4 different sizes.

Positioned on the lumbar region in order to trigger the patient's natural raising scheme, its rigid foam padding at the back-part level reduces painful pressure and its belt allows transfers in a sitting position.



	REFERENCES
Jersey	SWYP3-PF-S
Jersey	SWYP3-PF-M
Jersey	SWYP3-PF-L
Jersey	SWYP3-PF-XL
Padded Buttocks Belt	+PAD-PF

S	A	81	L	A	107
	B	21		B	29
M	A	94	XL	A	123
	B	25		B	33



Option:
padded buttocks belt

STAND-UP LIFTS SLINGS