




FAST ADJUSTABLE SLING

Comfort Model (without headrest)

The Fast Adjustable Sling's (without headrest) foam leggings reduce the pressures exercised on the thighs and the hooks with 2 settings (specifically designed for 2-point spreader bars) allow fast combinations for sitting, semi-sitting and lying transfers.

We recommend using it for bed to chair to bed transfers requiring a management of the head support (with the addition of removable headrest).



S	A	59
	B	62
	C	60
M	A	67
	B	68
	C	64
L	A	65
	B	78
	C	71
XL	A	70
	B	72
	C	90



REFERENCES	
Size S	SST-PF-S
Size M	SST-PF-M
Size L	SST-PF-L
Size XL	SST-PF-XL
Removable Headrest	+TAMOV-PF



Jersey




FAST ADJUSTABLE SLING

Bath Model (without headrest)

The Fast Adjustable Sling's (without headrest) foam leggings reduce the pressures exercised on the thighs and the hooks with 2 settings (specifically designed for 2-point spreader bars) allow fast combinations for sitting, semi-sitting and lying transfers.

We recommend using for transfers in a bath, a shower or a pool requiring a management of the head support (with the addition of removable headrest) as it is made in net fabric in order to evacuate the water and quickly dry.



S	A	58
	B	85
	C	65
M	A	67
	B	93
	C	76
L	A	74
	B	93
	C	83
XL	A	80
	B	100
	C	92



REFERENCES	
Size S	SSTBA-PF-S
Size M	SSTBA-PF-M
Size L	SSTBA-PF-L
Size XL	SSTBA-PF-XL
Removable Headrest	+TAMOV-PF



Net

