

# FAST ADJUSTABLE SLING

## Comfort Model




The Fast Adjustable Sling's integrated and boned headrest protects the cervical for patients without head support, its foam leggings reduce the pressures exercised on the thighs and the hooks with 2 settings (specifically designed for 2-point spreader bars) allow fast combinations for sitting, semi-sitting and lying transfers.


We recommend using it for bed to chair to bed transfers requiring a management of the head support.



<b>S</b>	A	58
	B	85
	C	65
<b>M</b>	A	67
	B	93
	C	76
<b>L</b>	A	74
	B	93
	C	83
<b>XL</b>	A	80
	B	100
	C	92



	REFERENCES
Size S	SML-PF-S
Size M	SML-PF-M
Size L	SML-PF-L
Size XL	SML-PF-XL



Jersey

# FAST ADJUSTABLE SLING

## Bath Model




The Fast Adjustable Sling's integrated and boned headrest protects the cervical for patients without head support, its foam leggings reduce the pressures exercised on the thighs and the hooks with 2 settings (specifically designed for 2-point spreader bars) allow fast combinations for sitting, semi-sitting and lying transfers.

We recommend using for transfers, requiring a management of the head support, in a bath, a shower or a pool as it is made in net fabric in order to evacuate the water and quickly dry.



<b>S</b>	A	58
	B	85
	C	65
<b>M</b>	A	67
	B	93
	C	76
<b>L</b>	A	74
	B	93
	C	83
<b>XL</b>	A	80
	B	100
	C	92



	REFERENCES
Size S	SMLBA-PF-S
Size M	SMLBA-PF-M
Size L	SMLBA-PF-L
Size XL	SMLBA-PF-XL



Net