




Jersey



Presented model:  
with sewed headrest



# HAMMOCK WITHOUT CROSSING SLING

## Comfort Model

The Hammock Without Crossing Sling is a hammock sling with a U-shape which avoids both hips tightness and the placement of the sling under the patient's buttocks, making the installation and removal comfortable for the patient and accommodating for the caregiver.

We recommend using it for bed to chair to bed transfers for amputees and for "lifting from the floor" transfers in a lying position (page 49).




<b>S</b>	A	41
	B	56
	C	61
<b>M</b>	A	48
	B	56
	C	68
<b>L</b>	A	51
	B	63
	C	76
<b>XL</b>	A	54
	B	65
	C	96

REFERENCES	
Size S	SUCPEHA-PF-S
Size M	SUCPEHA-PF-M
Size L	SUCPEHA-PF-L
Size XL	SUCPEHA-PF-XL
Removable Headrest	+TAMOV-PF
Sewed Headrest	+TAT-PF



Net



# HAMMOCK WITHOUT CROSSING SLING

## Bath Model

The Hammock Without Crossing Sling is a hammock sling with a U-shape which avoids both hips tightness and the placement of the sling under the patient's buttocks, making the installation and removal comfortable for the patient and accommodating for the caregiver.

We recommend using it for transferring amputees in a bath, a shower or a pool as it is made in net fabric in order to evacuate the water and quickly dry.



<b>S</b>	A	41
	B	56
	C	61
<b>M</b>	A	48
	B	56
	C	68
<b>L</b>	A	51
	B	63
	C	76
<b>XL</b>	A	54
	B	65
	C	96

REFERENCES	
Size S	SUCBAHA-PF-S
Size M	SUCBAHA-PF-M
Size L	SUCBAHA-PF-L
Size XL	SUCBAHA-PF-XL
Removable Headrest	+TAMOV-PF
Sewed Headrest	+TAT-PF