




# HAMMOCK WITHOUT CROSSING

## sling (without headrest)

The Hammock Without Crossing Sling is a hammock sling with a U-shape which avoids both hips tightness and the placement of the sling under the patient's buttocks, making the installation and removal comfortable for the patient and accommodating for the caregiver.

### COMFORT MODEL

We recommend using it for bed to chair to bed transfers for amputees and for "lifting from the floor" transfers in a lying position (with the addition of a headrest).



D3

Jersey



SCAN ME to learn more



Presented model: headrest version




#### REFERENCES

S: SUCPEHA-PF-S L: SUCPEHA-PF-L  
M: SUCPEHA-PF-M XL: SUCPEHA-PF-XL  
REMOVABLE HEADREST: +TAMOV-PF

### BATH MODEL

We recommend using it for transferring amputees in a bath, a shower or a pool as it is made in net fabric in order to evacuate the water and quickly dry.



D3

Net



SCAN ME to learn more



#### REFERENCES

S: SUCBAHA-PF-S L: SUCBAHA-PF-L  
M: SUCBAHA-PF-M XL: SUCBAHA-PF-XL  
REMOVABLE HEADREST: +TAMOV-PF