

HOLST ARM SLING



It is designed to prevent or to completely reduce the glenohumeral sub-luxation (source of pains) while allowing a daily setting up by the patients themselves, without changing the initial adjustment.

The elasticated strap exerts a vertical traction on the humerus by applying the force of the traction to the elbow and to the homolateral shoulder.



REFERENCES	
Left Side	EC-HSTG-PF
Right Side	EC-HSTD-PF

ARM SLINGS



Available in 2 models (with or without a dorsal strap), they are used in case of fracture, tendonitis or sprain to promote healing of the upper limb or shoulder.

Esthetical and comfortable, they are also equipped with a thumb rest.

The addition of a dorsal strap keeps the sling against the user's body (dorsal strap model only).



ARM SLING



ARM SLING WITH DORSAL STRAP

REFERENCES	
Size 1	ECB-01-PF
Size 2	ECB-02-PF
Size 3	ECB-03-PF

REFERENCES		
Arm sling with dorsal strap	Size 1	ECBSM1-PF
	Size 2	ECBSM2-PF
	Size 3	ECBSM3-PF

PATIENT POSITIONING