


# HAMMOCK WITHOUT CROSSING SLING


Comfort Model



Jersey




Presented model:  
with sewed headrest





# HAMMOCK WITHOUT CROSSING SLING

Bath Model



Net



PATIENT LIFTS SLINGS

The Hammock Without Crossing Sling is a hammock sling with a U-shape which avoids both hips tightness and the placement of the sling under the patient's buttocks, making the installation and removal comfortable for the patient and accommodating for the caregiver.

We recommend using the comfort model for bed to chair to bed transfers for amputees and for "lifting from the floor" transfers in a lying position (with the addition of a headrest).

We recommend using the bath model for transferring amputees in a bath, a shower or a pool as it is made in net fabric in order to evacuate the water and quickly dry.

REFERENCES		Comfort Model	Bath Model	Removable Headrest	+TAMOV-PF
	Size S	SUCPEHA-PF-S	SUCBAHA-PF-S	Sewed Headrest	+TAT-PF
Size M	SUCPEHA-PF-M	SUCBAHA-PF-M			
Size L	SUCPEHA-PF-L	SUCBAHA-PF-L			
Size XL	SUCPEHA-PF-XL	SUCBAHA-PF-XL			