

INTEGRAL BIRDY BELT



It is designed to prevent forward sliding of the patient and to maintain the bust in an upright position.

It can be installed while the patient is sitting to avoid any pelvic constriction.



We recommend using it if the patient's body has no holding on the armchair.

	REFERENCES
Size S	CBRIN2HP-PF-S
Size M	CBRIN2HP-PF-M
Size L	CBRIN2HP-PF-L
Size XL	CBRIN2HP-PF-XL

	A
Size S	60-80 cm
Size M	75-95 cm
Size L	90-110 cm
Size XL	105-130 cm

Maximum useful area: 200 cm



Y & PELVIC SLIM BELT



It is designed to prevent forward sliding of the patient and to maintain the bust in an upright position.

Its shoulder's straps can be separately adjusted in order to position the patient and to compensate from a bad position.



We recommend using it if the patient's body has no holding on the armchair.

	REFERENCES
Size S	CSYP2HP-PF-80
Size M	CSYP2HP-PF-95
Size L	CSYP2HP-PF-110
Size XL	CSYP2HP-PF-130

	A
Size S	60-80 cm
Size M	75-95 cm
Size L	90-110 cm
Size XL	110-130 cm

Maximum useful area: 230 cm

