



INTEGRAL BODY HOLDER

It is designed to prevent forward sliding of the patient and to maintain the bust in an upright position.

It is manufactured in a lamination of foam to provide a complete comfort.

We recommend using it if the patient's body has no holding on the armchair.



	A
Size S	60-80 cm
Size M	75-95 cm
Size L	90-110 cm
Size XL	110-130 cm



Maximum useful area: 230 cm

	REFERENCES
Size S	CBIM2HP-PF-S
Size M	CBIM2HP-PF-M
Size L	CBIM2HP-PF-L
Size XL	CBIM2HP-PF-XL



1. ABDOMINAL YZA BELT



3. PELVIC YZA BELT



2. BUST YZA BELT



4. INTEGRAL YZA BELT



YZA BELT

The YZA Belt is designed to provide the patient a customized chair support while avoiding the patient compression phenomena and is declined in 4 modular versions:

1. The Abdominal YZA Belt
2. The Bust YZA Belt is composed of the Abdominal YZA Belt to which is added a bust modular element
3. The Pelvic YZA Belt is composed of the Abdominal YZA Belt to which is added a pelvic modular element
4. The Integral YZA Belt is composed of the Abdominal YZA Belt to which is added a bust modular element and a pelvic modular element

We recommend using it if the patient's body has no holding on the armchair.



	REFERENCES
Abdominal YZA Belt	YZA2CAHP-PF
Pelvic YZA Belt	YZA2CAHP-PF+YZA2-MP-PF
Bust YZA Belt	YZA2CAHP-PF+YZA2-MB-PF
Integral YZA Belt	YZA2CAHP-PF+YZA2-MB-PF+YZA2-MP-PF

Maximum useful area: 200 cm