




LIFTING BELT



The Lifting Belt is used to raise a sitting person and preserves the caregiver of awkward postures during manual transfers.

Its inside fabric is in anti-slipping PU coating.



	REFERENCES
Size M	CREL-PF-M
Size L	CREL-PF-L



WALKING AID & LIFTING BELT



The Walking Aid & Lifting Belt is designed to prevent poor handling during manual transfers and is used on the patient as a walking aid for a safe handling.

As a lifting device, it is used around the patient's waist to provide a gripping point to the caregiver.

It is equipped with a self-centering 50 mm buckle to adjust the belt around the patient's waist.



	REFERENCES
Size M	CADM2-PF-S
Size L	CADM2-PF-M
Size M	CADM2-PF-L
Size L	CADM2-PF-XL