


# FLYER XL

With Electrical Opening Legs



## KEY FEATURES



Ergonomically designed 4-point spreader bar for an optimized weight distribution



10 000 N actuator for an increased maximum weight capacity (300 kg)



Strong castors with double ball bearing

## PRESENTATION

With its tremendous lifting power (10 000 N actuator) and its large electrical opening legs (up to 108 cm) which provides safety and comfort, the FLYER XL is a specialized model designed to improve the daily lives of bariatric patients (up to 300 kg) whose size and weight are challenging.

## OPTIONS



Lithium Battery for more efficiency:  
 - up to 84 cycles  
 - 2 years warranty



A mountable NAUSICAA weighing system to determine the patients weight (page 55)

## REFERENCES

With standard battery	FLYER-XL-E-PF
With lithium battery	FLYER-XL-E-LI-PF

