

PATIENT LIFTS SLINGS RANGE PRESENTATION

NAUSICAA Medical, number 1 French manufacturer of hoists slings since 1993, provides a complete range of stand-up lifts slings for all patients, in compliance with REACH legislation and conform to European Standard NF EN ISO 10535:2007.

For each of our patient lifts, we offer a suitable sling, compatible with all patient lifts on the market with a 2-point or 4-point spreader bar, to ensure that patient's mobility and rehabilitation are assisted in the most appropriate manner possible.

Our slings have been developed in partnership with health professionals in order to ensure the best possible support for the patient as well as a quick and easy installation for the caregiver.

CE CERTIFICATION

All our slings are CE certified in compliance with EU requirements arising to the product's manufacturer.

FABRIC

Our slings are manufactured in polyester either in jersey, a very easy to handle fabric which provides an optimized and comfortable positioning, or in net which is ideal for bath and toilet, avoids maceration and is easy to dry.

WASHING

All our slings are washable at 70°C.

REACH COMPLIANCE

Our slings follow the REACH regulation which aims to progressively remove the most dangerous chemical substances in the European Union.

OPTIONS

The addition of a removable or sewed or removable and boned headrest helps to take in charge the patient's head and protects the cervical.



PATIENT LIFTS SLINGS RANGE PRESENTATION

Size Chart

Weight / Sizes	1m50	1m55	1m60	1m65	1m70	1m75	1m80	1m85	1m90	+ 1m90
40 kg	S	S	S	S	S	S	M	M	M	M
45 kg	S	S	S	S	S	S	M	M	M	M
50 kg	S	S	S	S	S	S	M	M	M	M
55 kg	S	S	S	S	S	S	M	M	M	M
60 kg	M	M	M	S	S	S	M	M	M	M
65 kg	M	M	M	S	S	S	M	M	M	M
70 kg	M	M	M	M	M	M	M	M	M	M
75 kg	L	L	M	M	M	M	M	M	M	M
80 kg	L	L	L	L	M	M	M	M	M	M
85 kg	L	L	L	L	M	M	M	M	M	M
90 kg	L	L	L	L	L	L	M	M	M	M
95 kg	XL	XL	XL	XL	L	L	M	M	M	M
100 kg	XL	XL	XL	XL	L	L	L	L	L	L
110 kg	XL	XL	XL	XL	L	L	L	L	L	L
120 kg	XL	XL	XL	XL	XL	XL	L	L	L	L
130 kg	XL	XL	XL	XL	XL	XL	XL	XL	XL	XL
160 kg	2 XL	2 XL	2 XL	2 XL	2 XL	2 XL	XL	XL	XL	XL
190 kg	3 XL	3 XL	3 XL	3 XL	2 XL	2 XL	XL	XL	XL	XL
220 kg	3 XL	3 XL	3 XL	3 XL	3 XL	2 XL	2 XL	2 XL	XL	XL
250 kg	4 XL	4 XL	3 XL	3 XL	3 XL	3 XL	3 XL	2 XL	2 XL	2 XL
280 kg	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	3 XL	3 XL	3 XL	3 XL
320 kg	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL

A size colour-code of the edges has been put in place in order to quickly identify the size of each sling:

- to the size S corresponds the Yellow colour
- to the size M corresponds the Red colour
- to the size L corresponds the Green colour
- to the size XL corresponds the Blue colour
- to the size TU (one size fits) corresponds the Black colour

In order to choose the right sling, the caregiver shall consider the 3 following factors:

- Patient's height
- Patient's waist size
- Patient's thighs size

HEADREST



TAMOV



TABAL



TAMAX