



SCAN ME
to view all
devices

PATIENT LIFTS SLINGS

Presentation

Due to its high-quality standards in order to deliver the best product, NAUSICAA Medical is the number 1 French manufacturer of hoists slings since 1993.

We provide a complete range of patient lifts slings for all patients and conform to European Standard NF EN ISO 10535:2007.

Selecting the correct sling is paramount

Each patient is unique, and it is essential to select a sling that best suits the patient for comfort, safety and practicality.

Compatibility with all devices on the market

For each of our patient lifts, we offer a suitable sling, compatible with all patient lifts on the market with a 2-point or 4-point spreader bar, to ensure that patient's mobility and rehabilitation are assisted in the most appropriate manner possible.

CE Certification

All our slings are CE certified in compliance with EU requirements arising to the product's manufacturer.

Fabric

Our slings are manufactured in polyester either in jersey, a very easy to handle fabric which provides an optimized and comfortable positioning, or in net which is ideal for bath and toilet, as it avoids maceration and is easy to dry.

Washing

All our slings are washable at 70°C.

Reach Compliance

Our slings follow the REACH regulation which aims to progressively remove the most dangerous chemical substances in the European Union.



Certified



Quality

