



SIMPLE SLIM BELT

It is designed to prevent patients from falling as it is a support on the chair combined with a comfort belt.

This simple and robust system is common to all Slim belts.

We recommend using it when the patient needs to have the back pressed against the chair back and the buttocks at the bottom of the seat.



	A
Size S	60-80 cm
Size M	75-95 cm
Size L	90-110 cm
Size XL	110-130 cm



Maximum useful area: 230 cm

	REFERENCES
Size S	CCS2HP-PF-80
Size M	CCS2HP-PF-95
Size L	CCS2HP-PF-110
Size XL	CCS2HP-PF-130



LARGE ABDOMINAL BELT


It is designed to prevent patients from falling and to compensate slackening of abdominal muscles.

Its large area of contact with the patient prevents abdominal constriction.

We recommend using it when the patient needs to have the back pressed against the chair back and the buttocks at the bottom of the seat.



	A
Size S	60-90 cm
Size M	80-105 cm
Size L	95-125 cm
Size XL	115-145 cm



Maximum useful area: 230 cm

	REFERENCES
Size S	CCA2HP-PF-S
Size M	CCA2HP-PF-M
Size L	CCA2HP-PF-L
Size XL	CCA2HP-PF-XL