

STANDARD HAMMOCK SLING

Comfort Model

Made in France




The Standard Hammock Sling's indented cut on the sides makes it the preferred sling to let under the patient as it has no thickness which cause decubitus ulcers and its leggings comfort foam reduces the shearing effect.

We recommend using it for bed to chair to bed transfers where the sling remains on the chair.



	REFERENCES
Size S	SHAC-PF-S
Size M	SHAC-PF-M
Size L	SHAC-PF-L
Size XL	SHAC-PF-XL
Removable Headrest	+TAMOV-PF
Sewed Headrest	+TAT-PF

S	$\frac{A}{B}$	62 96
M	$\frac{A}{B}$	66 104
L	$\frac{A}{B}$	74 113
XL	$\frac{A}{B}$	83 121



Jersey

STANDARD HAMMOCK SLING

Bath Model

Made in France


The Standard Hammock Sling's indented cut on the sides makes it the preferred sling to let under the patient as it has no thickness which cause decubitus ulcers and its leggings comfort foam reduces the shearing effect.

We recommend using for transfers where the sling remains on the chair in a bath, a shower or a pool as it is made in net fabric in order to evacuate the water and quickly dry.



	REFERENCES
Size S	SHABA-PF-S
Size M	SHABA-PF-M
Size L	SHABA-PF-L
Size XL	SHABA-PF-XL
Removable Headrest	+TAMOV-PF
Sewed Headrest	+TAT-PF

S	$\frac{A}{B}$	62 96
M	$\frac{A}{B}$	66 104
L	$\frac{A}{B}$	74 113
XL	$\frac{A}{B}$	83 121



Net