




# STANDARD HAMMOCK

## sling

The Standard Hammock Sling's indented cut on the sides makes it the preferred sling to let under the patient as it has no thickness which cause decubitus ulcers and its leggings comfort foam reduces the shearing effect.

### COMFORT MODEL

We recommend using it for bed to chair transfers where the sling remains on the chair.



Jersey

**D2**



SCAN ME  
to learn more




#### REFERENCES

S: SHAC-PF-S                      L: SHAC-PF-L  
M: SHAC-PF-M                    XL: SHAC-PF-XL  
REMOVABLE HEADREST: +TAMOV-PF

### BATH MODEL

We recommend using for transfers where the sling remains on the chair in a bath, a shower or a pool as it is made in net fabric in order to evacuate the water and quickly dry.



Net

**D2**



SCAN ME  
to learn more



Presented model:  
with sewed headrest version

#### REFERENCES

S: SHABA-PF-S                      L: SHABA-PF-L  
M: SHABA-PF-M                    XL: SHABA-PF-XL  
SEWED HEADREST: +TAMOV-PF