



SCAN ME
to learn more

UNIVERSAL ECO

sling (comfort or bath)


The Universal Eco Sling's U-shape avoids the placement of the sling under the patient's buttocks, its leggings comfort foam reduces the shearing effect and the positioning of hooks at the shoulders level avoids the phenomenon of compression during a transfer.

The bath model is preferred for the bath or shower as it is made in net fabric in order to evacuate the water and quickly dry.

We recommend using it for bed to chair to bed transfers.



Jersey



Net

A1

REFERENCES

COMFORT-JERSEY: **SURPC-PF-TU**
BATH-NET: **SURP-PF-TU**
SIZE: **ONE SIZE FEET**
REMOVABLE HEADREST: **+TAMOV-PF**




SCAN ME
to learn more

UNIVERSAL ECO TOILET

sling

The Universal Eco Toilet Sling's extended U-shape completely frees the patient's buttocks to facilitate the positioning during the toilet and improve the installation and removal of the sling.

We recommend using it for toilet and daily change transfers and we do not advise it for patients with very low tone.



A1

Jersey



REFERENCES

JERSEY: **STORP-PF-TU**
BONED HEADREST: **+TABAL-PF**