



FLYER RANGE PRESENTATION



DESIGNED TO PERFORM HIGH QUALITY CARE

Practical and polyvalent, the FLYER range is our bestselling hoist line for several years as it includes high-end products for caregivers in order to provide quality care during transfers for all patients.

The FLYER range is declined in 4 different models:

- Standard sized with mechanical opening (page 44)
- Standard sized with electrical opening (page 45)
- Compact sized with mechanical opening (page 46)
- XL bariatric version with electrical opening (page 47)

The devices are suited for high-priority situations where fast and effective lifting is required.



LITHIUM BATTERY

As options, a lithium battery is available to enhance the use of the device and a mountable NAUSICAA weighing system (page 55) to determine the patient's weight.

NAUSICAA medical FLYER RANGE

SCAN
to watch



FLYER

With Mechanical Opening Legs



FLYER PLUS

With Mechanical or Electrical Opening Legs



FLYER COMPACT

With Mechanical Opening Legs



FLYER XL

With Electrical Opening Legs

FLYER RANGE KEY FEATURES

WARRANTY
5
YEARS



Equipped with advanced functional key features.

The FLYER range is designed to have the highest level of quality and safety to ensure that the caregiver and the patient are protected against any risk of injury.

With their exceptional very low base, these devices are compatible with the lowest medical beds on the market.

Versatile and useful, they are perfectly adapted to transfer a patient from a bed to a chair or a wheelchair.




01 4-POINT WIDE SPREADER BAR for an optimized weight distribution for patient comfort

Better for:

- Large or tall patients
- Patient positioning at 90° for transfer to a chair
- Includes head protection
- Lifting from floor



02 REMOTE AND CONTROL BOX
- 2 FUNCTIONS or 4 FUNCTIONS (when electrical opening legs is chosen) with battery level display on the remote control



03 OPENING LEGS:
- Mechanical
- Electrical



04 CASTORS:
- Reduced frame height (85 mm)
- 4 Double castors 75 mm
- 2 Single braked castors 100 mm

