

2 OR 4-POINT STANDARD SPREADER BAR



4-POINT WIDE SPREADER BAR



ADVANTAGES	<ul style="list-style-type: none"> - Small patients - Enveloping effect - Quick sling installation - Patient positioning at 90° to transfer to a chair 	<ul style="list-style-type: none"> - Patient comfort - Large patients - Patient positioning at 90° to transfer to a chair - Head protection - Lifting from floor
DISADVANTAGES	<ul style="list-style-type: none"> - Large patients - Patient comfort - Lifting from floor 	/

SPREADER BARS

Our patient lifts are equipped with two different spreader bars shapes depending on the devices.

Both models provide different advantages of use which are listed beside.

	LEAD BATTERY	LITHIUM BATTERY
Autonomy on a 35 cm lifting stroke actuator (FLYER - BLUE FLYER)	40 cycles	84 cycles
Time to recharge up to 85%	2h30	1h10
Lifespan (number of complete charges)	Up to 300	Minimum of 2000
Lifespan (number of micro charges)	Up to 240	Up to more than 8000
Weight of the battery pack	4 kg	1,8 kg
Warranty	None	2 years
Average lifetime	6 months to 2 years	2 years to 8 years

LITHIUM BATTERY OPTION

Lithium batteries support very well partial charging which reflects more the reality of use as you can recharge them as soon as the device is no longer used, even for a few minutes.

More efficient at all levels, and also lighter in handling, lithium batteries are perfectly suited for intensive use conditions.

