

FLYER

With Mechanical Opening Legs

Made in France
CE
 WARRANTY 5 YEARS

570/1800 mm
 85 mm
 660/1090 mm
 160 kg
 40 kg



KEY FEATURES



4-point spreader bar for an optimized weight distribution



6 000 N actuator for a 160 kg weight capacity



Reduced height of the front castors to pass under the lowest beds

PRESENTATION

Our best-seller provides a great ease of movement to perform secure and comfortable transfers for all patients (160 kg weight capacity). Its mechanical opening base legs (up to 109 cm) allows an ease of access to large widths wheelchairs.

OPTIONS



Lithium Battery for more efficiency:
 - up to 84 cycles
 - 2 years warranty



A mountable NAUSICAA weighing system to determine the patients weight (page 55)

REFERENCES

With standard battery	FLYER-PF
With lithium battery	FLYER-LI-PF