


TECHNICAL DATA FLYER RANGE



FLYER With Mechanical Opening Legs




DIMENSIONS

- A. Total length: 119 cm
- B. Minimum lifting height: 57 cm
- C. Maximum lifting height: 180 cm
- D. Total height: 132 cm
- E. Frame height: 8,5 cm
- F. Minimum base width: 66 cm
- G. Maximum base width: 109 cm
- H. Turning diameter: 133 cm



FLYER PLUS With Mechanical or Electrical Opening Legs




DIMENSIONS

- A. Total length: 119 cm
- B. Minimum lifting height: 57 cm
- C. Maximum lifting height: 180 cm
- D. Total height: 132 cm
- E. Frame height: 8,5 cm
- F. Minimum base width: 66 cm
- G. Maximum base width:
 - Mechanical Opening Legs: 109 cm
 - Electrical Opening Legs: 114 cm
- H. Turning diameter: 133 cm



FLYER COMPACT With Mechanical Opening Legs




DIMENSIONS

- A. Total length: 100 cm
- B. Minimum lifting height: 58 cm
- C. Maximum lifting height: 161 cm
- D. Total height: 120 cm
- E. Frame height: 8,5 cm
- F. Minimum base width: 52 cm
- G. Maximum base width: 89 cm
- H. Turning diameter: 110 cm



FLYER XL With Electrical Opening Legs



DIMENSIONS

- A. Total length: 125 cm
- B. Minimum lifting height: 55 cm
- C. Maximum lifting height: 182 cm
- D. Total height: 125 cm
- E. Frame height: 11 cm
- F. Minimum base width: 73 cm
- G. Maximum base width: 105 cm
- H. Turning diameter: 143 cm