

PATIENT LIFTS RANGE

Comparative



Weight Capacity
(in kg)



Total Weight
(in kg)



Minimum and
Maximum Lifting
Height (in cm)



Frame Height
(in cm)



Minimum and
Maximum Base
Width (in cm)

FLYER
Mechanical Opening Legs
page 38

FLYER COMPACT
Mechanical Opening Legs
page 39

FLYER PLUS
Mechanical or
Electrical Opening Legs
page 40

FLYER XL
Electrical Opening Legs
page 41

NAUSI-FLY
Mechanical Opening Legs
page 43

NAUSI-FLY COMPACT
Mechanical Opening Legs
page 44

NAUSI-FLY PLUS
Electrical Opening Legs
page 45

BLUE FLYER
Mechanical Opening Legs
page 47

BLUE FLYER COMPACT
Mechanical Opening Legs
page 48

BLUE FLYER PLUS
Electrical Opening Legs
page 49

	Weight Capacity (in kg)	Total Weight (in kg)	Minimum and Maximum Lifting Height (in cm)	Frame Height (in cm)	Minimum and Maximum Base Width (in cm)
FLYER Mechanical Opening Legs page 38	160	40	57 / 176	8	65 / 109
FLYER COMPACT Mechanical Opening Legs page 39	160	38	58 / 160	8	52 / 99
FLYER PLUS Mechanical or Electrical Opening Legs page 40	180	40 (mechanical) 42 (electrical)	59 / 180	8,5	65 / 109 (mechanical) 65 / 101 (electrical)
FLYER XL Electrical Opening Legs page 41	300	55	95 / 172	11	73 / 108
NAUSI-FLY Mechanical Opening Legs page 43	180	38	60 / 181	11,5	65 / 98
NAUSI-FLY COMPACT Mechanical Opening Legs page 44	160	31	61 / 160	11,5	57 / 88
NAUSI-FLY PLUS Electrical Opening Legs page 45	180	40	60 / 181	11,5	66 / 95
BLUE FLYER Mechanical Opening Legs page 47	180	38	60 / 181	11	65 / 98
BLUE FLYER COMPACT Mechanical Opening Legs page 48	160	31	61 / 160	11	57 / 88
BLUE FLYER PLUS Electrical Opening Legs page 49	180	40	60 / 181	11,5	66 / 95