

NAUSICAA WEIGHING SYSTEM

All our patient lifts can be equipped with a weighing system which is designed to easily determine the weight of patients up to 300 kg during transfers by a patient lift.

The device's maximum capacity on which it will be assembled prevails on the weighing system (delivered with fixations to reduce the loss of useful height).

It can only be mounted only on a NAUSICAA Medical patient lift and is supplied with fixings.

In order to maintain a correct useful height, we do not recommend to use it with compact devices such as the BLUE FLYER COMPACT, the NAUSI-FLY COMPACT and the MINI FLY'EVER.

FUNCTIONS:

TARE: turns the system to zero, even with a weight suspended to the system

ZERO: used to define the absolute zero

HOLD: displays a steady weight even if the system oscillates

UNIT: selects the desired unit of measure



SCAN ME
to learn more

REFERENCE
PESEE-NA-PF



300 kg



300 g



WARRANTY
1
YEAR