

U-SHAPED SLING

Comfort Model




The U-Shaped Sling's avoids the placement of the sling under the patient's buttocks, its leggings comfort foam reduces the shearing effect and the positioning of hooks at the shoulders level avoids the phenomenon of compression during a transfer.


We recommend using it for bed to chair to bed transfers and for "lifting from the floor" transfers in a lying position (page 49).



S	A	58
	B	54
	C	75
M	A	63
	B	58
	C	81
L	A	70
	B	61
	C	91
XL	A	77
	B	66
	C	102



	REFERENCES
Size S	SUCBAS-PF-S
Size M	SUCBAS-PF-M
Size L	SUCBAS-PF-L
Size XL	SUCBAS-PF-XL
Removable Headrest	+TAMOV-PF
Sewed Headrest	+TAT-PF



Jersey

U-SHAPED SLING

Bath Model




The U-Shaped Sling's avoids the placement of the sling under the patient's buttocks, its leggings comfort foam reduces the shearing effect and the positioning of hooks at the shoulders level avoids the phenomenon of compression during a transfer.

We recommend using it for transfers in a bath, a shower or a pool as it is made in net fabric in order to evacuate the water and quickly dry.



S	A	58
	B	54
	C	75
M	A	63
	B	58
	C	81
L	A	70
	B	61
	C	91
XL	A	77
	B	66
	C	102



	REFERENCES
Size S	SUCBAI-PF-S
Size M	SUCBAI-PF-M
Size L	SUCBAI-PF-L
Size XL	SUCBAI-PF-XL
Removable Headrest	+TAMOV-PF
Sewed Headrest	+TAT-PF



Net