

WAYUP RANGE PRESENTATION

RECOMMENDED FOR SEMI-ACTIVE TRANSFER:

The patient can lay both feet flat and has good bust support in order to be transferred in a motorized system while being actively stimulated.

The WAYUP range is declined in 3 different models (fixed base, mechanical or electrical opening) and is number 1 in sales for several years as it is an exceptional line by the combination of its compact design and the efficiency of its specific movement in one quick motion.

During a transfer, the devices allow the patient to reach a standing position by activating the bust tonicity while respecting the natural lifting scheme and supporting the lack of tonicity in the lower body.

As options, a lithium battery (page 12) is available to enhance the use of the device as well as a calf strap fixed to the shin block to provide a better feeling of safety.



WAYUP
fixed base
page 21



WAYUP
mechanical opening legs
page 22



WAYUP PLUS
electrical opening legs
page 23

OPTIONS



Lithium Battery



Calf Strap