

150 kg
37 kg
670/1600 mm
80 mm
340/560 mm



SCAN ME
to learn more

EASYLEV

fixed base

It has been perfected for patients with very low tonicity in the chest area and is equipped with a base adapted to a standard toilet shape, alongside its compactness with a 56 cm base width.

WATCH

the EASYLEV range product video
on our YOUTUBE channel



Traction belt for a patented combined movement increasing patient comfort



Premium control box for an enhanced use



Foam shin block for a comfortable support



Compact base adapted to toilet shape

OPTIONS



Lithium Battery



Calf Strap

REFERENCES

with standard battery: EASYLEV-2-PF
with lithium battery: EASYLEV-2-LI-PF