

FLYER COMPACT

With Mechanical Opening Legs

Made in France
CE
 WARRANTY 5 YEARS

580/1610 mm
 85 mm
 520/890 mm
 160 kg
 38 kg



KEY FEATURES



4-point spreader bar for an optimized weight distribution



Reduced height of the front castors to pass under the lowest beds



Easy transport with its quick folding capability

PRESENTATION

It is specially designed for home care use with its compact base and reduced dimensions, along with its mechanical opening legs (up to 99 cm) and its 160 kg weight capacity (6 000 N actuator).

OPTION



Lithium Battery for more efficiency:
 - up to 84 cycles
 - 2 years warranty

REFERENCES

With standard battery	FLYERCOMPACT-PF
With lithium battery	FLYERCOMPACT-LI-PF