



SCAN ME to learn more

# BACKREST sling

The Backrest Sling's U-shape associated with the net portion beneath the buttocks make the installation and removal comfortable for the patient and accommodating for the caregiver.

Its back part associated with the head-rest in rigid foam and the full stiffener allow to protect the spine during transfers.

We recommend using it for bed to chair to bed transfers and for "lifting from the floor" transfers in a lying position for specific pathologies.



A4

Jersey



### REFERENCES

S: SGLDOS-PF-S L: SGLDOS-PF-L  
M: SGLDOS-PF-M XL: SGLDOS-PF-XL




SCAN ME to learn more

# XXL slings

The following models, U-Shaped Sling - comfort and bath models (page 60), Quick Toilet Sling (page 61) and Hammock Without Crossing Sling - comfort model (page 65), are all available in several XL sizes (2 XL, 3 XL or 4 XL).

We recommend using them for bariatric patients transfers.



Jersey

U-Shaped Sling comfort model



A2

U-Shaped Sling bath model



A2

Quick Toilet Sling



C3

Hammock Without Crossing



D3

### REFERENCES

U-SHAPED SLING-JERSEY: SUCBAS-PF-2XL-3XL-4XL  
U-SHAPED SLING-BATH: SUCBAIT-PF-2XL-3XL-4XL  
BACKREST TOILET SLING: SUCPEHAT-PF-2XL-3XL-4XL  
QUICK TOILET SLING: STO-PF-2XL-3XL-4XL  
REMOVABLE HEADREST: +TAMAX-PF