



STAND-UP LIFTS SLINGS

Presentation



SCAN ME
to view all
devices

An historical know-how

Widely recognized as a textile expert, NAUSICAA Medical has been manufacturing slings for more than 25 years and works closely with health professionals to warrant that our slings perfectly match all patient needs.

Easy for the caregiver and safe for the patient

Quick and easy to install for the caregiver, our stand-up lifts slings are compliant to European Standard NF EN ISO 10535:2007 and ensure comfort as well as security for the patient.

One sling for each device

For each of our lifters, we offer an appropriate sling, compatible with all stand-up lifts on the market equipped with a hook system, to ensure that patient's mobility is assisted in the most appropriate manner possible for the safest transfer.

CE Certification

All our slings are CE certified in compliance with EU requirements.

Fabric

The slings are manufactured in polyester fabric.

Washing

All our slings are washable at 70°C.

Reach Compliance

Our slings follow the REACH regulation which aims to progressively remove the most dangerous chemical substances in the European Union.



Certified



Quality

