




BACKREST SLING

The Backrest Sling's U-shape associated with the net portion beneath the buttocks make the installation and removal comfortable for the patient and accommodating for the caregiver.

Its back part associated with the headrest in rigid foam and the full stiffener allow to protect the spine during transfers.

We recommend using it for bed to chair to bed transfers and for "lifting from the floor" transfers in a lying position for specific pathologies (page 49).



S	A	58
	B	90
	C	66
M	A	65
	B	96
	C	76
L	A	69
	B	99
	C	80
XL	A	74
	B	104
	C	89



REFERENCES	
Size S	SGLDOS-PF-S
Size M	SGLDOS-PF-M
Size L	SGLDOS-PF-L
Size XL	SGLDOS-PF-XL



XXL SLINGS

The following models, U-Shaped Sling - comfort and bath models (page 66), Quick Toilet Sling (page 67) and Hammock Without Crossing Sling - comfort model (page 71), are all available in several XL sizes (2 XL, 3 XL or 4 XL).

We recommend using them for transfers of bariatric patients.



U-Shaped Sling Comfort Model



U-Shaped Sling Bath Model



Hammock Without Crossing



Quick Toilet Sling



REFERENCES	
U-Shaped Sling-Jersey	SUCBAS-PF-2XL-3XL-4XL
U-Shaped Sling-Bath	SUCBAIT-PF-2XL-3XL-4XL
Hammock Without Crossing	SUCPEHAT-PF-2XL-3XL-4XL
Quick Toilet Sling	STO-PF-2XL-3XL-4XL
Removable Headrest	+TAMAX-PF

PATIENT LIFTS SLINGS