

# ENVELOPING SLING




The Enveloping Sling's integrated and boned headrest in addition to its enveloping shape make it suitable for the transfers of agitated and anxious patients.

Its U-shape avoids the placement of the sling under the patient's buttocks making the positioning more comfortable for the patient.

We recommend using it for bed to chair to bed transfers and for "lifting from the floor" transfers in a lying position (page 49).



<b>S</b>	A	60
	B	85
	C	97
<b>M</b>	A	73
	B	90
	C	110
<b>L</b>	A	74
	B	94
	C	124
<b>XL</b>	A	77
	B	102
	C	136



	REFERENCES
Size S	SLK-PF-S
Size M	SLK-PF-M
Size L	SLK-PF-L
Size XL	SLK-PF-XL



# ERGONOMIC SLING




The Ergonomic Sling's back part associated with its headrest in rigid foam lined with a 3D fabric increase comfort and allow to protect the spine during transfers.

Its high level of comfort and its integrated rigid headrest make it suitable for the transfers of agitated and anxious patients.

We recommend using it for bed to chair to bed transfers and for "lifting from the floor" transfers in a lying position (page 49).



<b>S</b>	A	57
	B	90
	C	60
<b>M</b>	A	60
	B	93
	C	70
<b>L</b>	A	72
	B	96
	C	80
<b>XL</b>	A	75
	B	103
	C	90



	REFERENCES
Size S	SMB-PF-S
Size M	SMB-PF-M
Size L	SMB-PF-L
Size XL	SMB-PF-XL



## PATIENT LIFTS SLINGS