




Jersey



FAST ADJUSTABLE ECO SLING (without headrest)

The Fast-Adjustable Eco Sling's foam leggings reduce the pressures exercised on the thighs and the hooks with 2 settings (specifically designed for 2-point spreader bars) allow fast combinations for sitting, semi-sitting and lying transfers.

We recommend using it for bed to chair to bed transfers requiring a management of the head.



TU	A	65
	B	69
	C	76

REFERENCES	
Jersey	SSTRP-PF-TU
Removable Headrest	+TAMOV-PF




Net



UNIVERSAL ECO HAMMOCK SLING

The Universal Eco Hammock's Sling indented cut on the sides makes it the preferred sling to let under the patient as it has no thickness which causes decubitus ulcers and it is preferred for the bath or shower due to its manufacturing in net fabric.

We recommend using it for "lifting from the floor" transfers in a lying position with the setup done at 180° by practicing a lateralization (page 49).



TU	A	74
	B	121

REFERENCES	
Net	HA-RBA-PF-TU
Removable Headrest	+TAMOV-PF

