

# UNIVERSAL ECO SLING

Comfort or Bath



# UNIVERSAL ECO TOILET SLING



Jersey



Net



Jersey



The Universal Eco Sling's U-shape avoids the placement of the sling under the patient's buttocks, its legging's comfort foam reduces the shearing effect and the positioning of hooks at the shoulders level avoids the phenomenon of compression during a transfer.

The bath model is preferred for the bath or shower as it is made in net fabric in order to evacuate the water and quickly dry.

We recommend using it for bed to chair to bed transfers.

REFERENCES	Comfort - Jersey	SURPC-PF-TU
	Bath-Net	SURP-PF-TU
	Removable Headrest	+TAMOV-PF

The Universal Eco Toilet Sling's extended U-shape completely frees the patient's buttocks to facilitate the positioning during the toilet and improve the installation and removal of the sling.

We recommend using it for toilet and daily change transfers and we do not advise it for patients with very low tone.

REFERENCES	Jersey	STORP-PF-TU
	Boned Headrest	+TABAL-PF