

PATIENT LIFTS RANGE Comparative



Weight Capacity
(in kg)



Total Weight
(in kg)



Minimum and
Maximum Lifting Height
(in cm)



Frame Height
(in cm)



Minimum and
Maximum Base
Width (in cm)

FLYER
Mechanical Opening Legs
(page 44)

FLYER PLUS
Mechanical or
Electrical Opening Legs
(page 45)

FLYER COMPACT
Mechanical Opening Legs
(page 46)

FLYER XL
Electrical Opening Legs
(page 47)

NAUSI-FLY
Mechanical Opening Legs
(page 52)

NAUSI-FLY
Electrical Opening Legs
(page 52)

NAUSI-FLY COMPACT
Mechanical Opening Legs
(page 53)

BLUE FLYER
Mechanical Opening Legs
(page 58)

BLUE FLYER
Electrical Opening Legs
(page 58)

BLUE FLYER COMPACT
Mechanical Opening Legs
(page 59)

	Weight Capacity (in kg)	Total Weight (in kg)	Minimum and Maximum Lifting Height (in cm)	Frame Height (in cm)	Minimum and Maximum Base Width (in cm)
FLYER Mechanical Opening Legs (page 44)	160	40	57 / 180	8,5	66 / 109
FLYER PLUS Mechanical or Electrical Opening Legs (page 45)	180	40 (mechanical) 42 (electrical)	57 / 180	8,5	66 / 109 (mechanical) 66 / 114 (electrical)
FLYER COMPACT Mechanical Opening Legs (page 46)	160	38	58 / 161	8,5	52 / 89
FLYER XL Electrical Opening Legs (page 47)	300	55	55 / 182	11	73 / 105
NAUSI-FLY Mechanical Opening Legs (page 52)	180	38	60 / 175	11	68 / 98
NAUSI-FLY Electrical Opening Legs (page 52)	180	40	60 / 175	11	68 / 108
NAUSI-FLY COMPACT Mechanical Opening Legs (page 53)	160	31	58 / 160	11	57 / 88
BLUE FLYER Mechanical Opening Legs (page 58)	180	38	65 / 180	11	68 / 98
BLUE FLYER Electrical Opening Legs (page 58)	180	40	65 / 180	11	68 / 108
BLUE FLYER COMPACT Mechanical Opening Legs (page 59)	160	31	63 / 165	11	57 / 88