

BLUE FLYER COMPACT

With Mechanical Opening Legs



KEY FEATURES



Securing shape of the spreader bar for the patient's head



Large handlebar for ease of handling



Reduced height of the front castors to pass under the lowest beds

PRESENTATION

It is specially designed for home care use with its 160 kg weight capacity (6 000 N actuator), its compact base, and its mechanical opening legs (up to 88 cm) which allows a better access to wheelchairs.

OPTION



Lithium Battery for more efficiency:
 - up to 84 cycles
 - 2 years warranty

REFERENCES

With standard battery	FLYERBLUE-CP-PF
With lithium battery	FLYERBLUE-CP-LI-PF

