

# NAUSIFLOW 2 System



130 kg  
WARRANTY 3 YEARS  
CE  
Decubitus Level

## NAUSIFLOW 2 Pump

It operates in alternating mode up to 130 kg, is intended for the prevention of very high-risk levels and is used as a care support for stage 3 pressure ulcers.

It has a standard pressure setting and its high flow pump allows to manage a wide range of mattresses.



### KEY FEATURES

- 01** Standard pressure settings: ability to manually adjust the pressure to optimize patient comfort
- 02** Self-management of seating positions (less than 2 hours)  
Management of seating positions (over 2 hours)
- 03** Power On/Off button
- 04** Visual and audible alarm: low pressure

## NAUSIFLOW 100-512 Air Mattress

Its compartmented configuration and its adjusted size (17 x 85 x 200 cm - 13 cm air height - 4 cm foam under-mattress) fits any type of 90 cm width bed (box spring from 85 to 90 cm).

It is equipped with an alternating heel section (4 cells) and a static head section (3 cells) to increase patient comfort.



Independent cells in polyurethane on nylon frame

Transport mode

CPR

REFERENCES	
NAUSIFLOW 2 Pump	NA2-COMP-PM
NAUSIFLOW 100-512 Mattress	NA512-MAT85-PM
NAUSIFLOW 2 System	NA2512-SYS