

-  180 kg
-  42 kg
-  620/1800 mm
-  113 mm
-  700/1070 mm




SCAN ME
to learn more

BLUE FLYER PLUS

electrical opening legs

It allows enhanced handling for daily transfers due to its flexibility of use (electrical openings up to 107 cm), its 8 000 N actuator (180 kg weight capacity), our premium control box and optimized castors for a better mobility.

WATCH

the BLUE FLYER range product video on our YOUTUBE channel



YouTube
SCAN to watch




360-degree rotation of the spreader bar for a simplified positioning



Premium control box for an enhanced use



8 000 N actuator offering an excellent compromise between maximum capacity (180 kg) and lifting speed



Electrical opening legs for ease of use

OPTIONS



Lithium Battery



Weighing System

REFERENCES

with standard battery: FLYERBLUE+-E-PF
with lithium battery: FLYERBLUE+-E-LI-PF