

DEMONSTRATION VIDEO

YouTube  
SCAN  
to watch



# BLUE WAYUP

## Fixed Base

Extremely compact with its 59 cm base and easy to handle, it is made for confined spaces, allowing easy access for the caregiver to perform a change during the toilet.



### KEY FEATURES



Fixed handles for a tight grip



Foam shin block for a comfortable support



Ball bearing castors for a better maneuverability

### OPTIONS




Lithium Battery for more efficiency:  
- up to 312 cycles  
- 2 years warranty



Calf Strap to secure the legs

### TECHNICAL DATA

- A. Total length: 94 cm
- B. Minimum lifting height: 79 cm
- C. Maximum lifting height: 123 cm
- D. Total height: 95 cm
- E. Frame height: 8 cm
- F. Minimum base width: 29,5 cm
- G. Maximum base width: 59 cm
- H. Turning diameter: 90 cm



### REFERENCES

With standard battery	WAYUPBLUE-PF
With lithium battery	WAYUPBLUE-LI-PF

STAND-UP LIFTS