

## PRODUCT SHEET NAUSIFLOW 2-512-E AIR MATTRESS



### Specifications:

Dimensions: 17 x 85 x 200 cm

13 cm air height

4 cm foam under-mattress

Top separable cover for easy maintenance -

CRIB 5 classified: antibacterial and antifungal

Class 1 medical device

### Key features:

Top separable cover

Separable zipper

Independent cells in polyurethane on nylon frame

Quick fasteners for the disconnectable and removable cells

In option, caps for cells discharge

Transport mode

CPR (cardiopulmonary resuscitation): deflates quickly (less than 20 seconds with a patient) the portion of the mattress which supports the patient's bust