

# NAUSIFLOW 2S System



130 kg

WARRANTY  
**3**  
YEARS



Decubitus  
Level

## NAUSIFLOW 2S Pump

It operates in alternating/static modes up to 130 kg, is intended for the prevention of very high-risk levels and is used as a care support for stage 3 pressure ulcers.

Equipped with a high-speed pump, it features a static operating mode to ensure treatment in palliative care.



### KEY FEATURES

- 01** Standard pressure setting: ability to manually adjust the pressure to optimize patient comfort
- 02** Palliative care mode: activates the static mode
- 03** Visual and audible alarm: low pressure required indicator
- 04** Self-management of seating positions (less than 2 hours)  
Management of seating positions (over 2 hours)

## NAUSIFLOW 100-512 Air Mattress

Its compartmented configuration and its adjusted size (17 x 85 x 200 cm - 13 cm air height - 4 cm foam under-mattress) fits any type of 90 cm width bed (box spring from 85 to 90 cm).

It is equipped with an alternating heel section (4 cells) and a static head section (3 cells) to increase patient comfort.



Independent cells in polyurethane on nylon frame



Transport mode



CPR

### REFERENCES

NAUSIFLOW 2S Pump	NA2S-COMP-PM
NAUSIFLOW 100-512 Mattress	NA512-MAT85-PM
NAUSIFLOW 2S System	NA2S512-SYS